

# Gamma-Sleep: Checklist Session 3 (sleep, experimental)

## Preparation

- Participant arrives, is led to same lab room as night before  
—> *note time of arrival in lab book*
- Dim the room lights
- Give time to complete usual sleep routine, ask to let us know when they're ready

## EEG setup

[SAME AS IN NIGHT 1]



Ask participant: anything else needed (toilet, water...)? If so, please do it now

- Attach mask to Neurofax, quick check if working → put on
- When participant is ready, lead them to bed, turn lights off  
—> *When fully ready, start recording*

## Overnight

[SAME AS IN NIGHT 1]

- As soon as EEG is scored as anything but W or N1 for 4 consecutive epochs (2 min) → start stimulation fade-in
- If PSG shows movement / wakefulness, open video stream → if indeed W (not just microarousal), fade out, then fade in again when back to deeper sleep (anything but W or N1)
- If participant must be awakened to adjust electrodes / mask during stimulation: fade out

## Morning

*[SAME AS IN NIGHT 1]*

- REDCap app: same except for
  - My projects —> gamma\_sleep —> Collect data —> GSQS\_exp
- after Sleep Quality Scale, thank for participation & terminate experiment
- participant payment form